

Analytical Chemistry

Experiment 5: Determination of Impurities in Whiskey Using Gas Chromatography

Your write-up for the lab is due by **Wednesday, May 2**—note change from syllabus




Chromatograms must be collected by Sunday, April 29

Each lab team should record in its notebook both what you do and observe when preparing the solutions on April 12. You do **not** need to write down your procedure when using the gas chromatograph. However, be sure to write down your observations.

Reference: Rice, Gary W., *J. Chem. Ed.* **1987**, *64*, 1055-1056.

Introduction

The uninitiated think that a distilled alcoholic beverage such as whiskey is simply an intense source of ethanol. However, the unique flavor of different brands of whiskey arises from a complex mixture of trace impurities. In this experiment you will determine the mixing ratios (in ppm) of three impurities found in two different whiskey samples using gas chromatography. The three analytes you will determine are shown below:



The relative amounts of these substances can serve as a fingerprint for a given brand of whiskey. You will use as your solvent a 40% solution of ethanol in water—this mimics the whiskey sample matrix. You will compensate for variations in sample injection volume by using an internal standard, 1-butanol.

Procedure

A. General Comments

1. During the lab session on April 12, each lab team will make up solutions. You will also need to sign up for a two-hour time slot to make measurements on the Hewlett-Packard 5890A Gas Chromatograph. (I will demonstrate its use the following week.) You should finish your measurements by April 29. (We will check out of lab the next day.)
2. Each stock bottle of analyte or internal standard has a dedicated 100- μ L syringe and a beaker. Pour a small amount of liquid into the beaker, and pipet out of the beaker. Please do not cross-contaminate the syringes, and please do not pipet directly out of any stock bottle.
3. There will be a stock solution of the 40% ethanol-in-water solvent.

B. Solution Preparation

You should partially fill each of the volumetric flasks with solvent before you add aliquot(s) of solute. (For the standards, the solvent is the 40% ethanol-in-water

solution. For the samples, the solvent is the whiskey itself.) The presence of solvent minimizes evaporative loss of solute, which would otherwise be a substantial source of systematic error.

1. *Preparation of whiskey samples.* Partially fill a 50-mL volumetric flask with one of the whiskey samples (be sure to note which one). You can fill the flask approximately halfway to the mark—the exact volume does not matter. Use a 100- μ L syringe to pipet a 50- μ L aliquot of 1-butanol into the same flask. Then dilute to the mark with the same whiskey sample. Mix thoroughly (you know the drill). Repeat for the second whiskey sample.
2. *Preparation of standard solutions.*
 - a. Put some solvent mixture (the 40% ethanol-in-water solution) into a clean 100-mL volumetric flask. Pipet 100 μ L of 1-butanol into the same flask. Dilute to the mark with the solvent mixture and mix thoroughly.
 - b. Repeat step a for each of the three analytes, 2-methyl-1-butanol, ethyl acetate, and 2-methyl-1-propanol. Each solution should have 100 μ L of one of the analytes.
 - c. Finally, make a 100-mL solution containing a 100- μ L aliquot of each of the three analytes, plus a 100- μ L aliquot of 1-butanol. That is, this solution should contain 400 μ L total of solute.

You should have a total of five standard solutions, each in a separate 100-mL volumetric flask.
3. For your data analysis, you will need to measure the density of the stock solution of ethanol in water. Perform a simple, yet accurate, experiment to determine this.

C. Using the HP 5890A Gas Chromatograph (GC)

Usually, the instructor or TA will have performed Steps 1-5. Skip these steps unless you are told otherwise.

1. Turn on the GC (the switch is on the back corner on the right hand side) and the integrator (the switch is on the lower right of the back panel.)
2. Open the three gas cylinders along the wall, using the main valves. The gas tank regulators are already set to the right pressures as long as you open the main valves completely. (As of Fall 2007, we use these pressures: H₂ at 30 psi, air at 30 psi, He at 40 psi.)
3. Start the flame ionization detector (FID) using the controls at the upper left of the instrument. The AUX GAS (helium) valve should be completely off. Open the AIR and HYDROGEN valves completely. Press the FID IGNITOR; you should hear a pop as the FID lights. The flame comes from the exothermic reaction of H₂(g) and O₂(g). To check that the FID is working, hold a watch glass over the detector vent; water vapor should condense on it. Then open the AUX GAS valve completely. Make sure the FID is still producing water vapor.

4. Turn on the integrator. Clean or replace, the ink cartridge. Hit the buttons ATT 2[^] (at the top of the keyboard) 0 ENTER (this sets the attenuation of the peaks to 2⁰ = 1).
5. Now program the GC. The injector and detector A port temperatures should be 200 °C. The oven should start at 60 °C for 5 minutes, then rise to 160° at a rate of 30°C/min. The helium gas (the mobile phase in the experiment) should have a flow rate of 100 mL/min. If this has not been set for you already, do the following: (a) Make sure COLUMN HEAD PRESSURE is set to 70 kPa. (b) Connect the soap film flow meter to the SPLIT VENT outlet. (c) Charge the flow meter with some fresh soap solution. (d) Squeeze the bulb at the bottom of the flow meter to get a bubble flowing. (e) Adjust the SPLIT VENT knob until a bubble rises to the 100-mL mark in 60 s.


Here is where your role begins:

- Start by obtaining a chromatogram of each of the four standard solutions containing one solute. The chromatogram should contain a huge peak for ethanol and a good-sized peak for the solute. (There may also be small peaks corresponding to trace impurities, and a big lump of unknown organic(s) eluting after ~10 min. Ignore these extra peaks.) This will let you determine the retention time of each solute.
- Then obtain a chromatogram of the standard containing all four solutes. (Based on your first three chromatograms, you will know how long you will have to wait to see all the peaks of interest.) This chromatogram will let you determine the response factor for each of the three analytes, as we will discuss in class.
- Finally, obtain chromatograms of your two whiskey samples.

Here are the specific procedures you should follow:

6. Sign the log book! Please note in the log book and in your lab notebook if you have any problems with the instrument.
7. Make sure the flame ionization detector (FID) is working by holding a watch glass over the detector vent. Water vapor should condense on it.
8. Since the GC has extremely low detection limits, it is important to minimize the cross contamination between runs. Therefore, before each run, flush the injection syringe 3 times with the solution you are about to analyze before injecting it onto the column. After each injection wash the syringe three times with the 40% ethanol-in-water solvent mixture. Use a Kimwipe to clean off any droplets sticking to the tip of the needle.
9. Use a small ($\leq 10\text{-}\mu\text{L}$) syringe to inject a 1- μL aliquot of the solution to be analyzed. (Never inject more than ~1- μL —it is easy to overload the narrow column.) Insert at least 1" of the needle into the injection port at the top left of the instrument. Simultaneously depress the syringe plunger and press the START button on the right side of the GC.
10. Press the TIME key on the control panel until it displays the ELAPSED time in minutes.

11. Let each elution continue until all expected substances have been detected. Then press the STOP button, and wait until the GC is ready before you start another run.
12. When you have completed all your runs, press the SHIFT and ENTER keys on the integrator simultaneously to advance the paper to the next perforation.
13. Keep the instrument on after you are done using it.



14. Note that on the chromatograms you will obtain (such as the one above), time is the vertical axis (increasing from top to bottom), and intensity is the horizontal axis (increasing from left to right). In the table at the bottom of the chromatogram, **RT** stands for retention time (in min), **AREA** gives the absolute peak area (which is proportional to the detector current), **WIDTH** gives the integrator's estimate of baseline peak width (in min) (so you can calculate resolution for two neighboring peaks if you wanted to), and **AREA%** gives the, well, percent area for each peak. (Ignore the **TYPE** column.) Use the **AREA** values as the peak intensities in your calculations of response factors and mixing ratios.

WASTE DISPOSAL: All solutions can go down the drain. As you are cleaning up, take the opportunity to make sure everyone in your lab team has every item on the cabinet inventory list. Each item should be clean and undamaged. This will facilitate lab check-out!

Data Analysis

The following calculations and discussion should all be presented in your notebook. **Also, please tape in all chromatograms into the notebook.** You may use a spreadsheet to facilitate your calculations, but this is not required. If you do use a spreadsheet, include a copy in your lab notebook and e-mail a copy to me.

1. Determine the retention time for each component by using the chromatograms you obtained for each of the four solutes by themselves. This will allow you to assign the peaks on the standard chromatogram containing all the analytes. Anticipate that the retention time for a given substance will not be completely precise.
2. Determine the response factor for each of the three analytes. Note that you will need to convert your experimental volume ratios to mass ratios (in ppm). This will require you to look up the density of each of the pure components and to use your experimental value for the solvent density. Be sure to write down the source of your density data.
3. Discuss the relationship (if any) between the analyte response factors you have measured and each analyte's molecular size and structure.
4. Determine the mixing ratios (in ppm) for 2-methyl-1-butanol, ethyl acetate, and 2-methyl-1-propanol in your two whiskey samples. You should assume that the mixing ratio of the 1-butanol internal standard in your two whiskey samples is equal to the 1-butanol mixing ratio in the standard solution. (This means that we are assuming that the density of all the solutions is the same.)
5. No quantitative error analysis is required for this experiment. However, you should provide a qualitative discussion of both systematic and random sources of error.